

Bed Rail Protector

Description:

Bed Rail Protectors are designed to protect patients from injuring themselves from hitting or getting limbs stuck between the bed rails, due to confusion or disorientation.

The Bed Rail Protectors are padded with foam for maximal patient safety and slip over the bed rail with two Velcro straps to fix them in place. They are made from PVC, with a removable cover that can be easily wiped clean with an antibacterial solution.



Dimensions:

Standard: 154cm × 48cm

Large: 172cm × 50cm

Instructions for use:

Slip the Bed Rail Protector over the bed rail and fasten in place with the two Velcro straps. Ensure the Bed Rail Protector is cleaned frequently.

The aid is a Class 1 Medical Device (non-invasive and non-sterile).