

Foot Plate Protectors

Description:

The Foot Plate Protectors provide a soft cushioning cover for hard wheelchair footplates. The covers fully encase the footplates to support and protect the patient, even when the footplates are raised. The Foot Plate Protectors will help to prevent pressure sores and increase patient comfort.

Hand or machine wash up to 80°C, and allow to dry naturally.

The Foot Plate Protectors come in a pair.

Dimensions:

Length: 23cm

Width: 18cm

Instructions for use:

The Foot Plate Protectors slide over the footplates and cover the hard surface entirely. The Foot Plate Protectors remains in place with the aid of the Velcro fastening strap.



The aid is a Class 1 Medical Device (non-invasive and non-sterile).