

Squeeze Handle Reacher



Ergonomic trigger handle

Clip attachment

360° rotation

Rubber lined jaws

The trigger mechanism can be pushed using the palm of the hand

Ideal for people with weakened fingers or hands, the BetterLiving Squeeze Handle Reacher helps people with mobility issues and assist in daily living tasks; helping live their life independently and in comfort.

The conventional trigger handle can be operated with the whole hand requiring less strength than the traditional two finger triggers; while the short trigger distance is ideal for stiff, arthritic fingers. The trigger also projects through the top of the handle, allowing the trigger mechanism to be pushed back using the palm of the hand.

For increased strength and durability the jaws are wire driven and the rubber lining provides a secure hold without damaging delicate items or fabric.

Rotating 360° the jaw allows for greater manoeuvrability and ease of access in hard to get places.

A magnet inserted into the tip of jaw assists with picking up small, hard to grab metal objects such as needles or hair pins. The reacher also has a removable clip allowing the user to easily attach the reacher to a walker or wheelchair .

Ideal for use in the home, garden or while out and about.

| | |
|--------------------|--------|
| CODE | BL0106 |
| LENGTH | 600 mm |
| JAW OPENING | 99 mm |
| WEIGHT | 200 g |
| COLOUR | Grey |
| WARRANTY | 1 year |

